


Leak Simulation Kit Instructions

The kit is easy to install and remove, creating a temporary “authentic” leak to test the efficacy of the LeakAlertor®.

Installation:

1. **Remove** the tank lid.
2. **Unwind** the weights from the strip
3. **Flush** toilet and let water drain from the tank. While holding the handle down, **Place** the small weight inside the drain approximately 2-3 inches. (See Figure 1)
4. **Release** the flush handle capturing the small weight inside the drain.
5. **Position** the large weight at the bottom of the tank approximately 2 inches from the drain.
6. **Bend** the strip over the edge of the tank to hang the leak simulator in place ensuring *light* tension on the string. (See Figure 2)
7. **Replace** tank lid.

The small weight is now partially down the drain with the string under the lip of the flapper. The string creates an obstruction preventing the flapper from seating properly, causing an authentic leak.



Removal:

1. **Remove** the tank lid.
2. **Grasp** the strip.
3. **Flush** the toilet, allow the water to drain, and then *carefully Lift* the simulator out of the drain.
4. **Replace** the tank lid.

While removing the weights, be careful not to allow them to break off from the string as they may cause a blockage or damage to the toilet.

We would appreciate hearing from you! Email comments to info@leakalserter.com. Please visit www.leakalserter.com for more information on our distribution programs.

For our limited warranty and more information, please visit us at: www.leakalserter.com or call us at 1-877-LEAKALERTOR (1-877- 532-5253)